FRASTEMA s.r.l.
ophthalmic equipments

**RIUNITO
NEW THREE STARS**

Art. 65 TA

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***NEW THREE STARS
UNIT***

INSTRUCTION MANUAL

COMPOSIZIONE RIUNITO NEW THREE STARS COMPOSITION OF NEW THREE STARS UNIT

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- 1 - Body of the unit
- 2 - Control panel
- 3 - Ophthalmoscope and skiascope housing
- 4 - Tables arm
- 5 - Instruments table
- 6 - Drawer cabinet
- 7 - Lamp columns and projector
- 8 - Standard chair

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
1 - BODY OF THE UNIT

The unit is composed of an internal welded structure (A) made with tubular iron where is welded the central support for the rotation of the instruments tables. Externally it's entirely covered by fiberglass repairs (B,C,D,E) that hosting the control panel and the housing for the instruments ophthalmoscope and skiascope.

Between repairs and frame is also inserted the drawer cabinet (F) with three drawers with one double.

For the transport of the unit is foreseen an hole in the anterior part of the repair to insert a tool (G) that go inside the frame. In the back part is foreseen a second lifting tool (H) to be inserted inside the columns holes.

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
The control panel is made in aluminium and lexan, is silk screen printed and has luminous indicators

- A: Main switch
- B: Projector control
- C: Auxiliary output socket control
- D: Halogen lamp control (*)
- E: Light room control
- F: Chair raising control
- G: Chair reset control
- H: Chair lowering control

(*) With a short touch switch on and switch off, keeping pushed adjust light intensity.

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(Fig. 2)

3 - OPHTHALMOSCOPE AND SKIASCOPE HOUSING

The ophthalmoscope and skiascope housing comprises two PVC containers with bushing for the connection of instruments, fed by low tension A.C. Current.

The output tension is controllable through the low tension output change over switch. (Fig. 1)

Upon request desk charger for ophthalmoscope and retinoscope.

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
Desk charger
for
ophthalmoscope
and retinoscope

Fig. 2

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
4 - TABLES ARM

The arm (A) is constructed in aluminium jointed to the shaft by means of needle bearing (C) fitted into cast iron support which allow the rotation (F).

There is a magnetic lock during operational use and a mechanical lock when the unit is out of operation. The arm built in a friction device for help operated to stop it in the work position (E). To make transport easier the arm with the table tops is fastened to the body of the unit using two set screw (B).

Adjusting the position of the arm (G) is possible to determinate the limit position of the tables during out of operation.

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5 - INSTRUMENTS TABLES

The structure is made in electrowelded mild steel plate with three wood table tops painted (1,2,3).


Table top nr. 1 is predisposed for the connection with the slit lamp. On it there are: jack socket (A), the rheostat for the adjustment of the light intensity (B) and the cables for the connection of the fixation point (6 V) (C).

Table top nr. 2 is predisposed for the connection of the ophthalmometer and on it there is only the jack socket.

Table top nr. 3 is predisposed for the connection of the autorefractometer or for other instruments powered at 220 V with a outlet type VDE (E).


Low tension instruments are automatically poered when are in working position. The outlet for the instrument powered at 220 V (110 V) is always in tension.

A magnetic system lock the arm and the table in everyone of the working position. The release is get pushing the push button (F) placed on the external side of every table top. On the same side is predisposed also the push button for the rising and the lowering of the chair (G).

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
6 - DRAWERS

The new three stars unit have a series of three services drawers which two with normal deep (A) and the other at double deep (B) where is possible to insert an instrument as tonometer.

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7 - LAMP AND PROJECTOR COLUMNS

Are composed with two steel columns (A and B) centered on two hollow pins (C) locked to the frame through them cross the electric wires. On the lamp column (A) is inserted a bended pipe (D) that go to the fiberglass lamp repair where is placed an halogen lamp (150W). On the same column is foreseen a series of holes where is possible to mount the counterbalanced arm (vertical type).

The projector column (B) has inserted a chromed pipe (F) adjustable. On the upper part is applied the projector

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8 - STANDARD CHAIR


It's consist of a shaped, curved support frame made of wood and has a covered anatomically designed seat and back support (A). The armrest are made of pressed polyurethane. The chair cover is in syntetic material in various colours. The chair can be rotated both to the right and to the left. As optional can be furnished with a tilting footrest (B).

The base (C) where lay is made with a shaped plate covered with a cover in preformed self-extinguishing polystyrene.

The hydraulic raising system (D) consist of a motor pump group (E) and an hydraulic press with the technical datas under indicated.

TECHNICAL FEATURES CHAIR RAISING HYDRAULIC SYSTEM

Power requirement	0.2 kW
Voltage	220V - 50Hz
Delivery	1.3-1.6 l/min
Thrust delivered	300 kg
Stroke	190 mm

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
OPERATING INSTRUCTIONS

1 - Connect the chair at the electric connector placed in the back of the unit

2 - Connect the power supply cable at the room electrical system. To be sure that this electrical system is protected by a ground fault interrupter.

3 - Switch on the main switch placed in the outlet then push the button POWER of the control panel.

Note : The reset position is adjustable. It's advisable to adjust the position sensor when the unit is installed, this adjustment allow to reduce the overall dimensions of the unit

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
Peso : kg 120
Potenza : 1110 VA
Tensione : 220-240 V (110 V)
Frequenza : 50-60 Hz

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
- 1- Chair connector
- 2- Auxiliary outlet 220 V (110V) 50/60Hz
- 3- Main switch

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- 1- Power supply transformer
- 2- Instruments table power supply transformer
- 3- On/off relay 12VAC 28.12
- 4- On/none/on relay 12VAC 28.13
- 5- Momentary relay 12VAC 5532
- 6- Lamp control switch and dimmer
- 7- Ophthalmoscope and retinoscope voltage commutator
- 8- Ophthalmometer voltage selector
- 9- Slit lamp voltage selector
- 10- PC board comp 04

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OPERATING OF THE UNIT

The New Three Stars unit is predisposed to utilize three ophthalmic instruments.

The position 1 is predisposed to utilize the sliding lamp alternating voltage powered (a commutator allow to select the voltage among 6-7.5-9-10.5-12-13.5 Volts). The switch on of the instrument happen automatically when the table is in working position. Two gray conductors powered at 6 V are used for the fixation point.

The standard preset voltage is of 7.5 V adjustable with the rheostat.

The position 2 is predisposed for the ophthalmometer powered with alternating voltage. A commutator allow to select the voltage among 6-7.5-9-10.5-12-13.5 Volts. The switch on of the instrument happen automatically when the table is in working position.

The standard preset voltage is of 6 V.

The position 3 is predisposed for the instrument powered at 220 V (autorefractometer, no contact tonometer, auto ref-keratometer etc.).

Powering of this instruments is interrupted from the main switch. On the external side of the table top support there is the button to unlock the positioning magnets and the switch for the rising and the lowering of the chair.

To place the choosen instrument in working position, rotate top tables on itself, stop matching the chair and approach until the lock position.

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MAINTENANCE

Equipment is maintenance free
External surfaces must be cleaned with detergent diluted in water
Do not use benzine, trichloroethylene, alcohol, solvents or acids.

NOTES FOR USES

Maximum admissible loading is 25 kg for every table top. To equilibrate the total load are used 8 counterweight weighting about 9 kg everyone. The columns are predisposed to utilize the phoropter arm both in traditional version than with vertical movement. The projector plate can be rotated and adjusted on the height

DESIGN USE

The unit New Three Stars art. 65/TA is a device designed to position one or more ophtalmic instruments at the correct position for visus examination

CLASSIFICATION

The unit New Three Stars art. 65/TA is an items of:

- Class I equipment
- Type B equipment
- Common equipment (without water proof casing)
- Equipment designed for continue operation
- Equipment not suitable for use in environments of flammable mixture of anaesthetics in air or oxigen or nitrous oxide